


















General Survey Ankylos C/X

Product line	X-Line		K-Line		Z-Line
Insert					
Distance height	4.40 mm	5.90 mm	2.95 mm	4.45 mm	
REF	I.03.04.X440	I.03.04.X590	I.03.04.K295	I.03.04.K445	
Denture magnet					
Distance height	2.65 mm		5.00 mm		
REF	U.00.01.X265R		U.00.01.K500		
Positioning cuff/ Resilience ring					
Resilience	0.3 mm	0.3 mm	0.00 mm	0.00 mm	0.3 mm
REF	P.00.01.X1	P.00.01.X2	P.00.01.K1	P.00.01.K2	P.00.05.K1
Impression post					
Distance height	6.95 mm		Please use the denture magnet for impression!		
REF	A.00.02.X695				
Laboratory replica		 Stainless steel			
Distance height	9.00 mm	9.00 mm	7.50 mm		
REF	M.00.01.X900	M.00.05.X900	M.00.01.K750		
Torque wrench insert	 <p>We offer a suitable torque wrench insert for the most popular ratchets.</p>				