








































General Survey XIVE/Frialit-2

Product line	X-Line	K-Line	Z-Line				
Insert ø 3.40 mm Distance height REF	 3.00 mm I.06.05.X300	 4.50 mm I.06.05.X450	 2.00 mm I.06.05.K200	 3.50 mm I.06.05.K350			
Insert ø 3.80 mm Distance height REF	 3.00 mm I.06.01.X300	 4.50 mm I.06.01.X450	 6.00 mm I.06.01.X600	 2.00 mm I.06.01.K200	 3.50 mm I.06.01.K350	 3.50 mm I.06.01.Z350	
Insert ø 4.50 mm Distance height REF	 3.00 mm I.06.02.X300	 4.50 mm I.06.02.X450	 6.00 mm I.06.02.X600	 2.00 mm I.06.02.K200	 3.50 mm I.06.02.K350	 3.50 mm I.06.02.Z350	
Insert ø 5.50 mm Distance height REF	 3.00 mm I.06.03.X300	 4.50 mm I.06.03.X450	 6.00 mm I.06.03.X600			 3.50 mm I.06.03.Z350	 5.00 mm I.06.03.Z500
Insert ø 6.50 mm Distance height REF						 3.50 mm I.06.04.Z350	 5.00 mm I.06.04.Z500

Denture magnet Distance height REF	 2.65 mm U.00.01.X265R	 5.00 mm U.00.01.K500	 3.15 mm U.00.01.Z315			
Positioning cuff/ Resilience ring Resilience REF	 0.3 mm P.00.01.X1	 0.3 mm P.00.01.X2	 0.00 mm P.00.01.K1	 0.00 mm P.00.01.K2	 0.3 mm P.00.05.K1	 0.3 mm P.00.03.Z1
Impression post Distance height REF	 6.95 mm A.00.02.X695	Please use the denture magnet for impression!	 6.95 mm A.00.02.Z695			
Laboratory replica Distance height REF	 9.00 mm M.00.01.X900	 9.00 mm M.00.05.X900 Stainless steel	 7.50 mm M.00.01.K750	 10.00 mm M.00.01.Z1000		
Torque wrench insert						
In addition to the torque wrench insert for Frialit 2 screwdriver 45/5722, we provide a torque wrench insert for the most popular ratchets.						