

































Titanmagnetics® inserts

Product line	X-Line ø 4.8 mm	Z-Line ø 5.8 mm	K-Line ø 5.2 mm	T-Line ø 5.8 mm
Surface				
	Fig. Height	Fig. Height	Fig. Height	Fig. Height
Ahead (BioComp Instruments)	 2.75 mm	 3.25 mm	 1.75 mm	 3.00 mm
External hexagon (Southern Implants IE, Brånemark, Ponto)	 3.50 mm 5.00 mm 6.50 mm 8.00 mm	 4.75 mm 6.00 mm 7.00 mm 8.00 mm	 3.00 mm 4.50 mm 6.00 mm	 5.20 mm 7.20 mm
EO (Straumann)	 3.25 mm 5.50 mm 7.25 mm 9.00 mm	 3.75 mm 6.50 mm 8.00 mm	 3.00 mm	 5.50 mm
Epitec-Grid direct (Stryker-Leibinger)	 3.00 mm 4.50 mm 6.00 mm 8.00 mm	 3.80 mm	 2.00 mm	
Epitec-Grid, Implant post (Stryker-Leibinger)	 4.20 mm			 4.50 mm
Ti-Epiplating, Protheses plates (Medicon)	 2.60 mm 4.50 mm 6.00 mm 8.00 mm	 4.50 mm	 2.50 mm 4.50 mm	 4.00 mm
Ti-Epiplating, Base post (Medicon)	 2.60 mm	 3.30 mm	 2.00 mm	 3.00 mm
Vistafix 3 (Cochlear)	 2.75 mm 4.50 mm	 4.00 mm 5.50 mm		 3.50 mm 5.00 mm
Torque wrench adapter	For all implant-systems fitting torque wrench inserts are available. For Z-Line inserts also torque wrench inserts of T-Line can be used.			



Titanmagnetics® Accessories

		Protheses magnets		Impression posts	Laboratory replicas	Healing Flanges
		universal	with collar			
X-Line ø 4.8 mm 	For	Acrylic	Acrylic		magnetic	
	Height/ring ø REF	2.65 mm U.00.01.X265R	2.65 mm U.00.01.X265K	6.95 mm A.00.02.X695	9.00 mm M.00.01.X900 16.00 mm M.00.01.X1600	4.20 mm/12.0 mm E.00.04.X1
Z-Line ø 5.8 mm 	For	Silicone	Silicone	Orbita	stainless steel	
	Height/ring ø REF	2.65 mm/11.5 mm U.00.02.X265	2.65 mm/11.5 mm U.00.02.X265K	5.55 mm A.00.01.X555	9.00 mm M.00.05.X900	
Z-Line ø 5.8 mm 	For	Acrylic	Acrylic			
	Height REF	3.15 mm U.00.01.Z315	3.15 mm U.00.01.Z315K	6.95 mm A.00.02.Z695	10.00 mm M.00.01.Z1000	
K-Line ø 5.2 mm 	For	Silicone	Silicone			
	Height/ring ø REF	3.15 mm/11.5 mm U.00.02.Z315	3.15 mm/11.5 mm U.00.02.Z315K			
K-Line ø 5.2 mm 	For	Acrylic				
	Height REF	5.00 mm U.00.01.K500		5.00 mm U.00.01.K500*	7.50 mm M.00.01.K750	
T-Line ø 5.8 mm 	For	Silicone				
	Height/ring ø REF	5.00 mm/11.5 mm U.00.02.K500				
T-Line ø 5.8 mm 	For	Acrylic				
	Height REF	5.70 mm U.00.01.T570		7.50 mm A.00.02.T750	10.50 mm M.00.01.T1050	
T-Line ø 5.8 mm 	For	Silicone				
	Height/ring ø REF	5.70 mm/11.5 mm U.00.02.T570				7.20 mm/12.0 mm E.00.04.T2

* For the impression use the facial prosthesis magnet.

