

With this I order acc. to the Steco General Terms and Conditions (see www.steco.de/en) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:
Please ask for shipping details

Customer-No.

Contact person

Date

Signature

Address (office or laboratory stamp):

CONTENT

Article	Planning	Pilot drill	Tube in tube	fully guided	Page
Titanium single sleeves	✓	✓	✗	✗	1
Titanium collar sleeve	✓	✓	✗	✗	1
Titanium reference balls	✓	✗	✗	✗	1
Titanium inner sleeves	✗	✓	✓	✗	2
Titanium inner sleeves with depth stopp	✗	✓	✓	✗	2
Titanium inner sleeve for endodontics	✗	✓	✓	✗	2
Titanium outer sleeves	✗	(✓)	✓	✗	2
Titanium outer sleeve open for lateral access	✗	(✓)	✓	✗	2
CeHa drill sleeves [inner]	✗	✓	✓	✗	3
CeHa drill sleeve [outer]	✗	(✓)	✓	✗	3
Titanium inner sleeves for Thommen Medical	✗	✓	✓	✗	3
Titanium outer sleeve for Thommen Medical	✗	(✓)	✓	✗	3
Titanium guided sleeves	✗	✗	✗	✓	4
Anchor pin sleeve					4

Fixation of mucosa-supported surgical guides.

SINGLE SLEEVES and Equipment

Titanium single sleeves (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
2.0	3.0	5.0	M.27.01.D200L5	10	83,00 €	
2.35	3.0	5.0	M.27.01.D235L5	10	83,00 €	
2.35	3.0	10.0	M.27.01.D235	10	83,00 €	

Accessories	Diameter	REF	PU	Price per PU	Quantity per PU
Template drill	3.0	M.27.01.B300	1	49,00 €	
Pressing tool	2.35	M.27.01.E235	1	19,50 €	

Titanium collar sleeve (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
2.0	3.0	5.0	M.27.31.D200L5	10	99,50 €	
2.2	3.0	5.0	M.27.31.D220L5	10	99,50 €	
2.35	3.0	5.0	M.27.31.D235L5	10	99,50 €	

REFERENCE BALLS

Titanium reference balls (All data in mm)

Diameter	REF	PU	Price per PU	Quantity per PU
5.0	M.27.09.D500	10	83,00 €	
2.5	M.27.09.D250	3	18,0 €	

Please pay attention to the different packaging units (PU)!

Discount already starts from an order size of 30 units (PU).

The different scales of discount are applied to the total order size of sleeves.

From 30 sleeves

10 % Discount

From 50 sleeves

15 % Discount

From 100 sleeves

25 % Discount



With this I order acc. to the Steco General Terms and Conditions (see www.steco.de/en) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:
Please ask for shipping details

Customer-No.

Contact person

Date

Signature

Address (office or laboratory stamp):

DOUBLE SLEEVES and Equipment

Titanium inner sleeves (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
1.5	3.5	6.0	M.27.03.D150L6	10	99,50 €	
1.5	3.5	10.0	M.27.03.D150	10	99,50 €	
1.6	3.5	6.0	M.27.03.D160L6	10	99,50 €	
2.0	3.5	6.0	M.27.03.D200L6	10	99,50 €	
2.0	3.5	10.0	M.27.03.D200	10	99,50 €	
2.1	3.5	6.0	M.27.03.D210L6	10	99,50 €	
2.2	3.5	6.0	M.27.03.D220L6	10	99,50 €	
2.25	3.5	6.0	M.27.03.D225L6	10	99,50 €	
2.35	3.5	6.0	M.27.03.D235L6	10	99,50 €	
2.35	3.5	10.0	M.27.03.D235	10	99,50 €	
2.5	3.5	6.0	M.27.03.D250L6	10	99,50 €	
2.8	3.5	6.0	M.27.03.D280L6	10	99,50 €	

Titanium inner sleeves with depth stop (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
1.16	3.5	5.0	M.27.24.D116L5	10	99,50 €	
1.3	3.5	5.0	M.27.24.D130L5	10	99,50 €	
2.0	3.5	5.0	M.27.24.D200L5	10	99,50 €	
2.2	3.5	5.0	M.27.24.D220L5	10	99,50 €	
2.3	3.5	5.0	M.27.24.D235L5	10	99,50 €	

Accessories	Diameter	REF	PU	Price per PU	Quantity per PU
Template drill	3.5	M.27.03.B350	1	73,00 €	
Pressing tool	2.0	M.27.03.E200	1	19,50 €	
Pressing tool	2.35	M.27.03.E235	1	19,50 €	

Titanium inner sleeve for endodontics (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
1.0	3.5	5.0	M.27.28.D100L5	10	99,50 €	

Titanium outer sleeves (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
3.5	4.0	5.0	M.27.02.D350L5	10	99,50 €	
3.5	4.0	6.0	M.27.02.D350	10	99,50 €	

Accessories	Diameter	REF	PU	Price per PU	Quantity per PU
Template drill	4.0	M.27.02.B400	1	73,00 €	
Pressing tool	3.5	M.27.02.E350	1	19,50 €	

Titanium outer sleeve open for lateral access (All data in mm)

Inner-ø	Outer-ø	Length	REF	VE	Price per PU	Quantity per PU
3.5	5.0	6.0	M.27.18.D350	5	79,50 €	

Please pay attention to the different packaging units (PU)!

Discount already starts from an order size of 30 units (PU).

The different scales of discount are applied to the total order size of sleeves.

From 30 sleeves

10 % Discount

From 50 sleeves

15 % Discount

From 100 sleeves

25 % Discount



With this I order acc. to the Steco General Terms and Conditions (see www.steco.de/en) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:
Please ask for shipping details

Customer-No.

Contact person

Date

Signature

Address (office or laboratory stamp):

DOUBLE SLEEVES (CeHa and Thommen)

CeHa drill sleeves [inner] (All data in mm)

Inner-ø	Outer-ø	Length	REF		PU	Price per PU	Quantity per PU
			Steco	Hafner			
1.6	4.5	5.0	M.27.06.D160	7814	10	102,50 €	
2.0	4.5	5.0	M.27.06.D200	7518	10	102,50 €	
2.1	4.5	5.0	M.27.06.D210	7930	10	102,50 €	
2.2	4.5	5.0	M.27.06.D220	7519	10	102,50 €	
2.35	4.5	5.0	M.27.06.D235	7520	10	102,50 €	
2.4	4.5	5.0	M.27.06.D240	7522	10	102,50 €	
2.5	4.5	5.0	M.27.06.D250	7523	10	102,50 €	
2.8	4.5	5.0	M.27.06.D280	7524	10	102,50 €	
3.0	4.5	5.0	M.27.06.D300	7525	10	102,50 €	
3.2	4.5	5.0	M.27.06.D320	7526	10	102,50 €	
3.3	4.5	5.0	M.27.06.D330	7527	10	102,50 €	
3.4	4.5	5.0	M.27.06.D340	7528	10	102,50 €	
3.5	4.5	5.0	M.27.06.D350	7529	10	102,50 €	
3.8	4.5	5.0	M.27.06.D380	7687	10	102,50 €	

Accessories	Diameter	REF		PU	Price per PU	Quantity per PU
		Steco	Hafner			
Template drill [inner]	4.5	M.27.06.B450	7533	1	73,00 €	

CeHa drill sleeves [outer] (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
4.5	5.2	5.0	M.27.05.D450	10	102,50 €	

Accessories	Diameter	REF		PU	Price per PU	Quantity per PU
		Steco	Hafner			
Template drill [outer]	5.2	M.27.05.B520	7534	1	73,00 €	

Titanium inner sleeves for Thommen Medical (VECTOdrill) (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
2.02	3.55	6.0	M.27.25.D200L6	10	99,50 €	
2.88	3.55	6.0	M.27.25.D280L6	10	99,50 €	

Titanium outer sleeve for Thommen Medical (VECTOdrill) (All data in mm)

Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
3.55	4.4	6.0	M.27.25.D350L6	5	63,75 €	

Please pay attention to the different packaging units (PU)!

Discount already starts from an order size of 30 units (PU).

The different scales of discount are applied to the total order size of sleeves.

From 30 sleeves

10 % Discount

From 50 sleeves

15 % Discount

From 100 sleeves

25 % Discount



With this I order acc. to the Steco General Terms and Conditions (see www.steco.de/en) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:
Please ask for shipping details

Customer-No.

Contact person

Date

Signature

Address (office or laboratory stamp):

GUIDED SLEEVES

Titanium guided sleeves (All data in mm)

For implant system/ surgery kit:	Diameter			SW	Length	REF	PU	Price per PU	Quantity per PU
	Inner	Outer	Collar						
Ankylos ND	4.48	5.5	6.5	5.1	4.0	M.27.15.D448	5	63,75 €	
Ankylos WD	4.9	5.6	6.6	5.2	4.0	M.27.15.D490	5	63,75 €	
Astra EV "ND"	4.6	5.6	6.6	5.2	4.0	M.27.15.D460	5	63,75 €	
Astra EV "WD"	5.2	6.0	7.0	5.6	4.0	M.27.15.D520	5	63,75 €	
Astra TX "small"	4.7	5.7	6.7	5.3	4.0	M.27.15.D470	5	63,75 €	
Astra TX "large"	5.7	6.7	7.7	6.3	4.0	M.27.15.D570	5	63,75 €	
BioHorizons Yellow	4.3	5.3	6.3	4.9	6.0	M.27.15.D432	5	63,75 €	
BioHorizons Green	5.1	6.0	7.0	5.6	6.0	M.27.15.D510	5	63,75 €	
BioHorizons Blue	6.3	7.0	8.0	6.6	6.0	M.27.15.D630	5	63,75 €	
Camlog ø 3.3	3.5	4.5	5.5	4.1	3.0	M.27.15.D350	5	63,75 €	
Camlog ø 3.8/4.3	4.5	5.5	6.5	5.1	3.0	M.27.15.D450L3	5	63,75 €	
Camlog ø 3.8/4.3	4.5	5.5	6.5	5.1	4.0	M.27.15.D450	5	63,75 €	
Camlog PL ø 3.3	3.6	4.7	5.5	5.2	3.0	M.27.15.D360L3	5	63,75 €	
Camlog PL ø 3.8/4.3	4.6	5.7	6.5	5.2	3.0	M.27.15.D460L3	5	63,75 €	
Camlog PL ø 5.0	5.6	6.7	7.5	6.2	3.0	M.27.15.D560L3	5	63,75 €	
MIS narrow	4.0	4.6	5.6	4.2	4.0	M.27.15.D400	5	63,75 €	
MIS	5.5	6.5	7.5	5.9	4.0	M.27.15.D550	5	63,75 €	
Nobel Biocare NP	4.1	4.7	5.7	4.3	3.5	M.27.15.D410	5	63,75 €	
Nobel Biocare RP	5.0	6.0	7.0	5.6	3.5	M.27.15.D500	5	63,75 €	
Nobel Biocare WP	6.2	7.0	8.0	6.6	3.5	M.27.15.D620	5	63,75 €	
PrimeTaper EV (XND)	3.4	5.4	4.0	4.0	4.0	M.27.15.D340	5	63,75 €	
PrimeTaper EV (ND)	4.6	5.6	6.6	5.2	4.0	M.27.15.D460	5	63,75 €	
PrimeTaper EV (WD)	5.2	6.0	7.0	5.6	4.0	M.27.15.D520	5	63,75 €	
Simplant universal RP	4.2	5.0	6.0	4.6	4.0	M.27.33.D420L4	5	63,75 €	
Simplant universal WP	5.2	6.0	7.0	5.6	4.0	M.27.33.D520L4	5	63,75 €	
Straumann Guided Pilot	2.2	3.8	4.8	3.4	6.0	M.27.26.D220L6	5	63,75 €	
Straumann	2.8	3.8	4.8	3.4	6.0	M.27.15.D280	5	63,75 €	
Straumann	5.0	6.0	7.0	5.6	5.0	M.27.15.D500L5	5	63,75 €	
Straumann iGuide	3.9	5.2	5.9	4.6	5.0	M.27.33.D390L5	5	63,75 €	
Straumann iGuide	5.2	6.5	7.2	5.9	5.0	M.27.33.D520L5	5	63,75 €	
Sweden & Martina GS415	4.25	5.0	5.5	4.6	5.0	M.27.15.D425	5	63,75 €	
Sweden & Martina GS550	5.55	6.35	7.0	5.95	5.0	M.27.15.D555	5	63,75 €	
OmniTaper EV (ND)*	4.48	5.5	6.5	5.1	4.0	M.27.15.D448	5	63,75 €	
OmniTaper EV (WD)*	5.2	6.0	7.0	5.6	4.0	M.27.15.D520	5	63,75 €	
Zimmer A	4.2	5.2	6.2	4.6	4.0	M.27.15.D420	5	63,75 €	
Zimmer B	5.3	6.3	7.3	5.7	4.0	M.27.15.D530	5	63,75 €	

*Formerly XIVE

SW = Wrench size

The implant names mentioned are trademarks of the respective manufacturers.

Anchor pin sleeve (All data in mm)

Implant system	Inner-ø	Outer-ø	Length	REF	PU	Price per PU	Quantity per PU
Nobel Biocare	1.5	3.5	10.0	M.27.20.D150L10	10	99,50 €	

Please pay attention to the different packaging units (PU)!

Discount already starts from an order size of 30 units (PU).

The different scales of discount are applied
to the total order size of sleeves.

From 30 sleeves

10 % Discount

From 50 sleeves

15 % Discount

From 100 sleeves

25 % Discount

